

# PATHOLOGY WORKSTATION AIP31-1800



# PATHOLOGY WORKSTATION AIP31-1800

Workstations are used during dissections, autopsies, necropsies, and other pathological procedures. It includes grossing/trimming stations and dissection tables. It may also include a built-in sink and ventilation system. Used in Hospital pathology department, Pathology laboratory.

Also known as Pathological sampling bench.

## AIP31-1800 PATHOLOGY WORKSTATION

Odor-remove function

Single water sink, two-person operations

Exhaust air system: High-power suction device, low noise

Optional water heater is available, realizing hot & cold water adjustment

The base cabinet is inward-concave type, convenient for sitting operation

Optional pulverizer can completely smash the waste during test, to avoid blocking drainage pipes

Automatic flushing device, slightly right-side inclined worktable, easy for waste discharge and cleaning



## SPECIFICATIONS

Model	AIP31-1800
Work Surface Height	850 mm
Exhaust Airflow Volume	500~600m <sup>3</sup> /h
Main Body Material	304 stainless steel
Work Table Material	304 stainless steel
Package Size (WxDxH)	1960x890x2170 mm
Gross Weight	340 kg
Working Area Size (WxDxH)	1695x430x900 mm
External Size (WxDxH)	1800x750x1950 mm
Consumption	300 W
Power Supply	AC220V±10%, 50/60Hz; 110V±10%, 60Hz



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)