

UV STERILIZATION CABINET AIP71-O253



UV STERILIZATION CABINET AIP71-O253

UV Sterilization Cabinet are used to sterile various equipments.They are manufactured to the highest standards to ensure reliable, sterile results. It rapidly disinfect clinical equipment, furniture and electronic devices.

Used in Hospital , Laboratory , House.

Also known as UV Disinfection Cabinet.

AIP71-O253 UV STERILIZATION CABINET



SPECIFICATIONS

Model	AIP71-O253
Type	With Ozone
UV Ray Wavelength	253.7/184.9nm
Sterilization Way	2 sides of UV ray+1 side of Ozone
Capacity	200 L
Power	120 W
UV Lamp Material	Quartz glass
Shelves	Standard 2pc, can be optional
Internal Size	518x523x778 mm
External Size	750x663x1200 mm
Packing Size	875x788x1265 mm
Gross Weight	100 kg



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com