

# UV STERILIZATION CABINET AIP71-W253



# UV STERILIZATION CABINET AIP71-W253

UV Sterilization Cabinet are used to sterile various equipments.They are manufactured to the highest standards to ensure reliable, sterile results. It rapidly disinfect clinical equipment, furniture and electronic devices.

Used in Hospital , Laboratory , House.

Also known as UV Disinfection Cabinet.

## AIP71-W253 UV STERILIZATION CABINET



### SPECIFICATIONS

Model	AIP71-W253
Type	Without Ozone
UV Ray Wavelength	253.7nm
Sterilization Way	3 sides of UV ray
Capacity	200 L
Power	120 W
UV Lamp Material	Quartz glass
Shelves	Standard 2pc, can be optional
Internal Size	518x523x778 mm
External Size	750x663x1200 mm
Packing Size	875x788x1265 mm
Gross Weight	100 kg



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)