

AUTO COAGULATION ANALYZER ANA13-160



AUTO COAGULATION ANALYZER ANA13-160

Auto Coagulation Analyzer used to test clotting items, Chromogenic, Immunologic measuring methods
Used in Clinical diagnostic, Laboratory, Pharmaceutical.
Also known as coagulation analyzer.

ANA13-160 AUTO COAGULATION ANALYZER



SPECIFICATIONS

Model	ANA13-160
Throughput	160T/H
Analysis Method	Coagulation method, Immunoturbidimetry
Analysis Items	PT, APTT, TT, FIB, IDD, FDP
Sample Positions	6
Tilt reagent Positions	11 refrigerated reagent positions



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@biozef.com | Website: biozef.com