

AUTOSAMPLER FOR TOC ANALYZER ANA55-20P



AUTOSAMPLER FOR TOC ANALYZER ANA55-20P

TOC analyzer directly shows the total organic carbon content in water or liquid solutions. It is designed to make testing easier, reproducible and accurate. TOC analyzers maximize both sensitivity and productivity, making them the ideal choice for monitoring in wide range of aqueous samples.

Used in Water for injection, Cleaning validation, Environmental, Chemical, Pharmaceutical Industries, Research, Laboratory.

Also known as Laboratory TOC Analyzer, Total Organic Carbon Analyzer.

ANA55-20P AUTOSAMPLER FOR TOC ANALYZER

Compact design, easy to put on bench top.

Easy to use, no need for special training.

20 positions auto sampling which saves time.

Adopts USB to connect with PC.

Connect with TOC to realize auto-sampling.

Equipped with Multi-sensors for real-time monitoring to ensure instrument safety.

Run for long time with good stability, no need for special maintenance.

Compatible with model AT121 and AT122.



SPECIFICATIONS

Model	ANA55-20P
Sample Position	20
Sample Volume	30/60 ml
Interface	USB
Output	120 W
Power Supply	100~240 V
Frequency	50/60 Hz



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com