

HORIZONTAL AUTOCLAVE AUT22-0185B



HORIZONTAL AUTOCLAVE AUT22-0185B

Laboratory horizontal autoclaves are useful in bulk sterilization applications of equipment and materials in various industries and research institutions. Also for disposing biological waste, reducing the risk of contamination and sterile the waste before discarding. We provide a variety of horizontal autoclave according to applications, features and customers requirements.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

AUT22-0185B HORIZONTAL AUTOCLAVE

Sterilizing programs for fabric, instrument, rubber and liquid

Touch screen display of temperature, pressure, time, running status and failure warning

Post-drying function makes the residual humidity of sterilized articles less than 3%

The positive pressure air-exhaust system to exhaust the cold air from chamber

Equipped with printer to print date, time and process parameters

Build-in efficient steam generator, microcomputer temperature control

Perfect four-level authority management system to avoid misoperation

Double door optional



SPECIFICATIONS

Model	AUT22-0185B
Capacity	185 L
Chamber Size	Φ500x950 mm
Chamber Material	SUS304
Max Working Pressure	0.23 Mpa
Working Temp.	115-138°C
Temperature Accuracy	+0.1° C
Consumption	11 KW
Power Supply	380 V, 50 HZ
External Size(LxWxH)	1350x750x1742 mm
Package Size(LxWxH)	1510x890x2060 mm
Gross Weight	540 kg



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com