

HORIZONTAL AUTOCLAVE AUT22-0241DA



HORIZONTAL AUTOCLAVE AUT22-0241DA

Laboratory horizontal autoclaves are useful in bulk sterilization applications of equipment and materials in various industries and research institutions. Also for disposing biological waste, reducing the risk of contamination and sterile the waste before discarding. We provide a variety of horizontal autoclave according to applications, features and customers requirements.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

AUT22-0241DA HORIZONTAL AUTOCLAVE



Product Image Coming Soon

SPECIFICATIONS

Model	AUT22-0241DA
Capacity	241 L
Chamber Size (Wx H x L)	600x600x740 mm
External Size (WxHxL)	370x1880x950 mm
Net Weight	580 Kg
Power Consumption	2.5 KW
Power Supply	380 V 50 Hz



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com