

HORIZONTAL AUTOCLAVE AUT22-0241SA




HORIZONTAL AUTOCLAVE AUT22-0241SA

Laboratory horizontal autoclaves are useful in bulk sterilization applications of equipment and materials in various industries and research institutions. Also for disposing biological waste, reducing the risk of contamination and sterile the waste before discarding. We provide a variety of horizontal autoclave according to applications, features and customers requirements.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

AUT22-0241SA HORIZONTAL AUTOCLAVE



Product Image Coming Soon

SPECIFICATIONS

Model	AUT22-0241SA
Capacity	241 L
Chamber Size (Wx H x L)	600x600x740 mm
External Size (WxHxL)	1370x1880x1000 mm
Net Weight	660 Kg
Power Consumption	2.5 KW
Power Supply	380 V 50 Hz



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com