

# HORIZONTAL AUTOCLAVE AUT22-0300B



# HORIZONTAL AUTOCLAVE AUT22-0300B

Laboratory horizontal autoclaves are useful in bulk sterilization applications of equipment and materials in various industries and research institutions. Also for disposing biological waste, reducing the risk of contamination and sterile the waste before discarding. We provide a variety of horizontal autoclave according to applications, features and customers requirements.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

## AUT22-0300B HORIZONTAL AUTOCLAVE

Sterilizing programs for fabric, instrument, rubber and liquid

Touch screen display of temperature, pressure, time, running status and failure warning

Post-drying function makes the residual humidity of sterilized articles less than 3%

The positive pressure air-exhaust system to exhaust the cold air from chamber

Equipped with printer to print date, time and process parameters

Build-in efficient steam generator, microcomputer temperature control

Perfect four-level authority management system to avoid misoperation

Double door optional



## SPECIFICATIONS

Model	AUT22-0300B
Capacity	300 L
Chamber Size	Φ632x1000 mm
Chamber Material	SUS304
Max Working Pressure	0.23 Mpa
Working Temp.	115-138°C
Temperature Accuracy	+0.1° C
Consumption	17 KW
Power Supply	380 V, 50 HZ
External Size(LxWxH)	1400x890x1780 mm
Package Size(LxWxH)	1700x1050x2090 mm
Gross Weight	830 kg



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)