

# HORIZONTAL AUTOCLAVE AUT22-0352DDB




# HORIZONTAL AUTOCLAVE AUT22-0352DDB

Laboratory horizontal autoclaves are useful in bulk sterilization applications of equipment and materials in various industries and research institutions. Also for disposing biological waste, reducing the risk of contamination and sterile the waste before discarding. We provide a variety of horizontal autoclave according to applications, features and customers requirements.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

## AUT22-0352DDB HORIZONTAL AUTOCLAVE



Product Image Coming Soon

### SPECIFICATIONS

Model	AUT22-0352DDB
Capacity	352 L
Chamber Size (Wx H x L)	600x600x980 mm
External Size (WxHxL)	1215x1880x1190 mm
Net Weight	1000 Kg
Power Consumption	24 KW
Power Supply	380 V 50 Hz



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)