

# HORIZONTAL AUTOCLAVE AUT22-0810DB



# HORIZONTAL AUTOCLAVE AUT22-0810DB

Laboratory horizontal autoclaves are useful in bulk sterilization applications of equipment and materials in various industries and research institutions. Also for disposing biological waste, reducing the risk of contamination and sterile the waste before discarding. We provide a variety of horizontal autoclave according to applications, features and customers requirements.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

## AUT22-0810DB HORIZONTAL AUTOCLAVE



Product Image Coming Soon

### SPECIFICATIONS

Model	AUT22-0810DB
Capacity	810 L
Chamber Size (Wx H x L)	610x910x1460 mm
External Size (WxHxL)	1310x2070x1650 mm
Net Weight	1260 Kg
Power Consumption	4 KW
Power Supply	380 V 50 Hz



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)