

HORIZONTAL AUTOCLAVE AUT22-0810SDB



HORIZONTAL AUTOCLAVE AUT22-0810SDB

Laboratory horizontal autoclaves are useful in bulk sterilization applications of equipment and materials in various industries and research institutions. Also for disposing biological waste, reducing the risk of contamination and sterile the waste before discarding. We provide a variety of horizontal autoclave according to applications, features and customers requirements.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

AUT22-0810SDB HORIZONTAL AUTOCLAVE



Product Image Coming Soon

SPECIFICATIONS

| | |
|-------------------------|-------------------|
| Model | AUT22-0810SDB |
| Capacity | 810 L |
| Chamber Size (Wx H x L) | 610x910x1460 mm |
| External Size (WxHxL) | 1310x2070x1700 mm |
| Net Weight | 1460 Kg |
| Power Consumption | 61 KW |
| Power Supply | 380 V 50 Hz |



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com