

# HORIZONTAL AUTOCLAVE AUT22-1500SB



# HORIZONTAL AUTOCLAVE AUT22-1500SB

Laboratory horizontal autoclaves are useful in bulk sterilization applications of equipment and materials in various industries and research institutions. Also for disposing biological waste, reducing the risk of contamination and sterile the waste before discarding. We provide a variety of horizontal autoclave according to applications, features and customers requirements.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

## AUT22-1500SB HORIZONTAL AUTOCLAVE



Product Image Coming Soon

### SPECIFICATIONS

Model	AUT22-1500SB
Capacity	1500 L
Chamber Size (Wx H x L)	680x1180x1870 mm
External Size (WxHxL)	1450 x1990x2130 mm
Net Weight	2211 Kg
Power Consumption	5 KW
Power Supply	380 V 50 Hz



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)