

# CLASS N MEDICAL AUTOCLAVE AUT35-03



# CLASS N MEDICAL AUTOCLAVE AUT35-03

Medical Autoclave assure proper cleaning and disinfection of instruments and utensils for medical, dental and laboratory field. It minimizes the risk of exposure to infection for each patient and provider. Due to its door locking system and small utilization of space, it becomes an ideal choice for the customers.

Used in Medical, Dental, Laboratory, Pharmaceutical, Research.

Also known as Medical Sterilizer, Dental Autoclave, Medical Autoclave, Pre and Post Vacuum Autoclave, Benchtop Autoclave, Benchtop Sterilizer, Benchtop Steam Sterilizer, Laboratory Medical Autoclave.

## AUT35-03 CLASS N MEDICAL AUTOCLAVE



Product Image Coming Soon

### SPECIFICATIONS

Model	AUT35-03
Capacity	3 L
Sterilization Pressure	0.22 Mp
Sterilization Temperature	134 °C
Sterilization Time	4 min
Sterilization Cycle time	11 min
Required water	Distilled water or purified water
Product Size	300X190X350 mm
Outer Dimension	370x260x460 mm
Gross Weight	7.5 Kg
Power	0.5 KW
Voltage	AC220/110V+10%



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)