

# DRY BATH INCUBATOR BAC1C-100RT1



# DRY BATH INCUBATOR BAC1C-100RT1

Compact design with stable operation, easy cleaning, replacement and disinfection of metal blocks. It is widely used for enzyme reactions, immunoassays, melting and boiling points, incubation and activation of cultures and variation of laboratory procedure. Flexible to accommodate a variety of interchangeable heating blocks.

Used in Calibration tests, High temperature organic distillation, External circulation, High temperature experiments. Also known as Bath Circulator, Laboratory Bath Circulator.

## BAC1C-100RT1 DRY BATH INCUBATOR

Microprocessor temperature controller

Easy maintenance

LCD display with touch key panel for temperature setting

Over temperature alarm system



### SPECIFICATIONS

Model	BAC1C-100RT1
Temperature Range	RT +5-100°C
Display Accuracy	0.1°C
Temperature Control Accuracy	±0.3°C

### OPTIONAL ACCESSORIES

Accessory Code	Name	Description
LS54689	Control panel (including board)	
LS54701	Heating tube	



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)