

# ANALYTICAL BALANCE BAL11-300



# ANALYTICAL BALANCE BAL11-300

Analytical Balances are designed to meet the specific weighing requirement in the laboratory working environment. Performance, durability, robust, easy-to-use, precise results, efficient handling, space saving makes analytical balance an ideal choice.

Used in Totalization, Formulation, SQC, Industry, Laboratory, Research.

Also known as Lab Balance, Laboratory Analytical Balance, Analytical Scale, Laboratory Analytical Scale.

## BAL11-300 ANALYTICAL BALANCE

Aluminum shell, Anti-jamming, Anti-radiation

The multiple functions modes with tare, counting, unit conversion and weighing

It is a magnetic sensor.



## SPECIFICATIONS

Model	BAL11-300
Capacity (g)	300 g
Resolution (g)	0.001 g
Min Weighing (g)	0.004 g
Stable Time	≤3 sec
Display	LCD(white back light with black font)
Sensor	Magnet
Pan Size	Φ 80 / Φ 90 mm
Operation Temperature	5-25°C
Repeatability	±0.002 g
Liner	±0.003 g
Cal. Weight	200 g
Overall Dimension	240x190x265 mm
Packing Size	Inside: 415x285x430 mm Outside: 475x315x485 mm
Gross Weight	10 kg

## OPTIONAL ACCESSORIES

Accessory Code	Name	Description
LS53461	RS232C & RS485 Interface	
LS53473	Printer	





**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)