

# ANALYTICAL BALANCE BAL12-120



# ANALYTICAL BALANCE BAL12-120

Analytical Balances are designed to meet the specific weighing requirement in the laboratory working environment. Performance, durability, robust, easy-to-use, precise results, efficient handling, space saving makes analytical balance an ideal choice.

Used in Totalization, Formulation, SQC, Industry, Laboratory, Research.

Also known as Lab Balance, Laboratory Analytical Balance, Analytical Scale, Laboratory Analytical Scale.

## BAL12-120 ANALYTICAL BALANCE

Electromagnetic balanced sensors

The multiple functions modes with tare, counting, unit conversion and weighing

External Calibration



## SPECIFICATIONS

Model	BAL12-120
Capacity (g)	120 g
Resolution (g)	0.0001 g
Min Weighing (g)	0.0004 g
Stable Time	≤3 sec
Display	LCD(white back light with black font)
Pan Size	Φ 80 mm
Operation Temperature	5-25°C
Repeatability	±0.0002 g
Liner	±0.0003 g
Cal. Weight	100 g
Overall Dimension	240x190x265 mm
Packing Size	Inside: 415x285x430 mm Outside: 475x315x485 mm
Gross Weight	10 kg

## OPTIONAL ACCESSORIES

Accessory Code	Name	Description
LS53461	RS232C & RS485 Interface	
LS53473	Printer	





**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)