

ANALYTICAL BALANCE BAL12-220



ANALYTICAL BALANCE BAL12-220

Analytical Balances are designed to meet the specific weighing requirement in the laboratory working environment. Performance, durability, robust, easy-to-use, precise results, efficient handling, space saving makes analytical balance an ideal choice.

Used in Totalization, Formulation, SQC, Industry, Laboratory, Research.

Also known as Lab Balance, Laboratory Analytical Balance, Analytical Scale, Laboratory Analytical Scale.

BAL12-220 ANALYTICAL BALANCE

Electromagnetic balanced sensors

The multiple functions modes with tare, counting, unit conversion and weighing

External Calibration



SPECIFICATIONS

Model	BAL12-220
Capacity (g)	220 g
Resolution (g)	0.0001 g
Min Weighing (g)	0.0004 g
Stable Time	≤3 sec
Display	LCD(white back light with black font)
Pan Size	Φ 80 mm
Operation Temperature	5-25°C
Repeatability	±0.0002 g
Liner	±0.0003 g
Cal. Weight	200 g
Overall Dimension	240x190x265 mm
Packing Size	Inside: 415x285x430 mm Outside: 475x315x485 mm
Gross Weight	10 kg

OPTIONAL ACCESSORIES

Accessory Code	Name	Description
LS53461	RS232C & RS485 Interface	
LS53473	Printer	



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@biozef.com | Website: biozef.com