

PRECISION BALANCE BAL26-4100



PRECISION BALANCE BAL26-4100

Precision Balance incorporates unique high end technology which is built to match the highest international standards. Outstanding weighing performance and innovative sample handling technologies are used. Due to its efficiency, error-free, and compliant features it becomes an unique product.

Used in Chemistry, Biological Research, Industry, Laboratory.

Also known as Loading Balance, Laboratory Precision Balance, Laboratory Loading Balance, Precision Scale, Laboratory Precision Scale.

BAL26-4100 PRECISION BALANCE

Aluminium shell, Anti-jamming, Anti-radiation

LED digital display

The multiple functions modes with tare, counting, unit conversion and weighing



SPECIFICATIONS

Model	BAL26-4100
Capacity (g)	4100 g
Resolution (g)	0.01 g
Min Weighing (g)	0.04 g
Stable Time	3 sec
Display	LED / LCD
Sensor	Load Cell
Pan Size	180x180 mm
Operation Temperature	5-25°C
Repeatability	±0.03 g
Cal. Weight	2000 g
Gross Weight	5 kg

OPTIONAL ACCESSORIES

Accessory Code	Name	Description
LS53461	RS232C & RS485 Interface	
LS53473	Printer	
LS53485	USB	



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@biozef.com | Website: biozef.com