

WEIGHING SCALE BAL33-510



WEIGHING SCALE BAL33-510

Weighing Scale uses high precision weighing sensor to keep weighing stable. Used in laboratory, industrial and food retailing application. Due to its latest design and friendly material, it is widely used in the laboratory.

Used in Weighing, Percent Weighing, Industrial, Food Retailing, Laboratory, Research.

Also known as Weighing Balance, Laboratory Weighing Scale, Laboratory Weighing Balance.

BAL33-510 WEIGHING SCALE

Stainless steel Pan

LCD digital display

The multiple functions modes with tare, counting, unit conversion and weighing



SPECIFICATIONS

Model	BAL33-510
Capacity (g)	510 kg
Resolution (g)	50 g
Min Weighing (g)	250 g
Stable Time	≤2 sec
Display	LED
Sensor	Load Cell
Pan Size	600x800 mm
Operation Temperature	5-25°C
Repeatability	±100 g
Liner	±150 g

OPTIONAL ACCESSORIES

Accessory Code	Name	Description
LS53461	RS232C & RS485 Interface	
LS53473	Printer	



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@biozef.com | Website: biozef.com