

# TABLETOP CENTRIFUGE CEN11-05



# TABLETOP CENTRIFUGE CEN11-05

These centrifuges are widely used in the fields of biochemistry, pharmacy, clinical medicine and other research fields for variety of separation needs. The programmable microprocessor centrifuge allows modification of the parameter during operation. In addition, accommodates various type of tubes or plates and size of the sample.

Used in Cell Separation, Precipitation, Sample Processing, Clinical, Cell Culture, Microplate Processing, Biochemistry, medical diagnosis.

Also known as Floor Standing Centrifuge, Laboratory Floor Type Centrifuge, Benchtop Centrifuge, Non Refrigerated High Speed Centrifuge, Laboratory Tabletop Centrifuge.

## CEN11-05 TABLETOP CENTRIFUGE

Gentle braking at low speeds

LCD display

Short Spin key for quick centrifugation

Lid is opened automatically after the operation is finished, which avoids overheating of the sample and save processing time

Safety lid-lock



## SPECIFICATIONS

Model	CEN11-05
Maximum Capacity (No of tubes x Vol.)	6x15ml
Maximum Speed RPM	100-5000
Maximum RCF	3074xg
Timer Range	15 S~99 min
Acceleration/ Deceleration	20-80s↑ 20-80s↓
Power	40 W
Power Supply	110V 50/60Hz, 230V-50/60Hz

## OPTIONAL ACCESSORIES

Accessory Code	Name	Description	Capacity
LS55158	Rotor	D12x75 glass tube	6x15ml
LS55159	Rotor	Conical tube with screw cap	8x15 ml
LS55160	Rotor	Conical tube with screw cap	12x5 ml





**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)