

# PCR CENTRIFUGE CEN32-0080PCR



# PCR CENTRIFUGE CEN32-0080PCR

These laboratory multi purpose Medical Centrifuges are ideal solution for high-end research experiment such as PRP/PPP, cell washing, HCT, blood type test etc. in the biology, radio-immunity, medical and other research fields. These versatile centrifuges have precise temperature control with highly efficient double loop control design which allows fast heating and cooling.

Used in Molecular Biology Separations, Clinical, Labs, Biochemistry, Tissue Culture, Genomics, Drug discover. Also known as Benchtop Centrifuge, Laboratory Centrifuge.

## CEN32-0080PCR PCR CENTRIFUGE

Compact design with small footprint

Acceleration / Deceleration: <18 seconds

Low noise

Metal rotor

Short Spin key for quick centrifugation

Temperature sensitive samples are protected by excellent air-flow

Lid is opened automatically after the operation is finished, which avoids overheating of the sample and save processing time

Maintenance free, Brushless DC motor

Autoclavable rotor



## SPECIFICATIONS

|                            |                            |
|----------------------------|----------------------------|
| Model                      | CEN32-0080PCR              |
| Maximum Speed RPM          | 800-13, 400                |
| Maximum RCF                | 12100xg                    |
| Timer Range                | 15 S~99 min                |
| Acceleration/ Deceleration | 13s / 12 s                 |
| Overall Dimension          | 260x260x150 mm             |
| Weight (Net/Gross)         | 4.3 kg                     |
| Power                      | 100 W                      |
| Power Supply               | 110V 60Hz, 230V 50 / 60 Hz |

## ACCESSORIES

| Accessory Code | Name  | Description                    | Max. Speed |
|----------------|-------|--------------------------------|------------|
| LS55200        | Rotor | 8x1.5/2.0ml, 8x0.2ml, 2 strips | 4000 rpm   |

## OPTIONAL ACCESSORIES

| Accessory Code | Name    | Description | Capacity        |
|----------------|---------|-------------|-----------------|
| LS55156        | Adapter |             | 0.2 ml PCR tube |
| LS55157        | Adapter |             | 0.5 ml PCR tube |





**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)