

# MICROCENTRIFUGE CEN43-060



# MICROCENTRIFUGE CEN43-060

This compact, powerful micro centrifuge is useful to separate or sediment small volume samples at high speed. The short-spin key, soft brake, FastTemp features optimizes the overall function of the centrifuges. Faults are diagnosed automatically and displays fault code. These hybrid micro-centrifuge can run PCR tubes, Capillary tubes along with standard tubes.

Used in Molecular Biology Separations, Clinical, Labs, Biochemistry, Tissue Culture, Genomics, Drug discover.  
Also known as Benchtop Centrifuge, Laboratory Centrifuge.

## CEN43-060 MICROCENTRIFUGE

Quick release of rotor for easy rotor exchange

Gentle braking system for quick sample processing



## SPECIFICATIONS

Model	CEN43-060
Maximum Capacity (No of tubes x Vol.)	8x1.5 ml / 2.0 ml
Maximum Speed RPM	4000 & 6000
Maximum RCF	1200xg & 2000xg
Overall Dimension	135x147x120 mm
Weight (Net/Gross)	0.7 kg
Power Supply	AC100-240V

## OPTIONAL ACCESSORIES

Accessory Code	Name	Description	Capacity
LS55142	Fixed-Angle Rotor		8x1.5/2.0ml
LS55143	Strip Rotor		4-8x0.2ml
LS55144	Adapter		0.5-0.6ml tube
LS55145	Adapter		0.2ml tube





**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)