

# STANDARD FUME HOOD FUH14-1500



# STANDARD FUME HOOD FUH14-1500

Created to provide protection from toxic fumes and hazardous volatile vapours, and are reliable and easy to use. Inbuilt ventilation system lowers the risk of fire and explosion. For metal analysis, pharmacy and research field along with safety measures, this product becomes unique choice for the buyers to purchase according to their needs.

Used in Metal analysis, Pharmacy, Research Field, Commercial, Medical, Laboratory.

Also known as Standard Fume Hood, Laboratory Fume Hood, Chemical Fume Hood.

## FUH14-1500 STANDARD FUME HOOD

UV lamp for sterilization

Adjustable air speed: 9 levels

Motorized front window, height adjustable

Microprocessor control system, LED display

With memory function in case of power failure

Cold-rolled steel with anti-bacteria powder coating

High grade melamine board with acid and alkali resistance function

Chemical resistant phenolic resin



## SPECIFICATIONS

Model	FUH14-1500
Airflow Velocity	0.3~0.8 m/s
Air Tightness	No leakage
Max Opening	520 mm
Length	1500 mm
LED Lamp	16Wx1
UV Lamp	Emission of 253.7 nanometers for most efficient decontamination
Front Window	5 mm toughened glass, Motorized, Height Adjustable
Exhaust Duct	4 meters PVC exhaust duct, $\phi$ 300 mm
Pipe Strap	2
Blower Optional Choice	Built-in centrifugal blower, Speed adjustable
Work Surface Height	850 mm
Inner Dimension	1320x670x730 mm
Overall Dimension	1540x800x2200 mm
Package Size (Main Body)	1690x990x1690 mm
Package Size (Base Cabinet)	1690x970x980 mm
Noise Level	$\leq$ 60 dB (A)
Weight (Net/Gross)	350 kg
Power	500 W
Power Supply	110/220V $\pm$ 10%, 50/60Hz

## ACCESSORIES

Accessory Code	Name	Description
LS52354	water tap	
LS52367	Gas tap	
LS52380	Water sink	
LS52393	Base cabinet	

Accessory Code	Name	Description
LS52406	Pipe strap	
LS52419	Total load of two waterproof sockets: 500W	
LS53479	UV Lampx2	
LS54931	LED lamp	

## OPTIONAL ACCESSORIES

Accessory Code	Name	Description
LS52445	Active carbon filter	



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)