

AUTOMATIC AIR GENERATOR GEN13-1000



AUTOMATIC AIR GENERATOR GEN13-1000

Air Generator is designed for removing hydrocarbons from the air stream and produces pure, hydrocarbon free and dry air, with a very low dew point. Long service life, easy maintenance, abrasion resistant, cost effective, performance, quick and easy servicing are some of the features which makes it convenient for the laboratory.

Used in Pharmaceutical, Gas Chromatography Detectors, Petrochemical, Laboratory, Medical, Research, Industries, Environmental, Clinical, Forensics..

Also known as Laboratory Air Generator.

GEN13-1000 AUTOMATIC AIR GENERATOR

The hydrocarbon is changed into carbon dioxide and vapor in the catalyzer with Pt-Pd carrier

Increased gas output quality

Compact design, easy using, automatically operation, safety



SPECIFICATIONS

Model	GEN13-1000
Hydrogen purity	99.999%
Flow Rate	0-1000 ml / min
Output Pressure	0-0.4 Mpa
Pressure Stability	<0.001 Mpa
Dimension	460x370x360 mm
Net Weight	20 kg
Consumption Power	400 W
Power Supply	220 V \pm 10%, 50 Hz



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com