

GLASSWARE WASHER GLA11-320



GLASSWARE WASHER GLA11-320

Lab Automatic Glassware washer is an automatic machine for the process of washing, disinfection and drying of various laboratory equipment such as flasks, beakers, pipettes and so on. Washing carried out through air heater and HEPA filter facilitates by drying

Used in University , Research institution, Pharmaceutical, Chemical industry, Microbiology Lab.

Also known as Washer Disinfectant.

GLA11-320 GLASSWARE WASHER

High efficient cleaning system, designed with European reliable pump, optimized spray arm and nozzles and self-cleaning program.

Rapid and efficient drying system, designed with independent air heater, HEPA filter and drying process

Water heating temperature can reach 99°C , while hot air drying temperature can reach 120°C

Professional configuration design-stainless steel for corrosion resistance; reinforced glass window mounted to achieve a clear view on washing process

PLC touch screen control with present 12 programs and 99 customized programs

Safety protection-electronic security door lock to prevent unexpected door opening , water & drying air temperature dual control and equipped with emergency switch.



SPECIFICATIONS

| | |
|----------------------------------|--|
| Model | GLA11-320 |
| Capacity | 320 L |
| Water Consumption Cycle | 22 L |
| Inlet Water Pressure | 0.3~0.8 Mpa |
| Noise | ≤ 55 dB |
| Material | Internal Material: SS316, Outer Material:SS304 |
| Overall Power Consumption | 6/17 KW |
| Power Consumption- Water Washing | 1.5 KW |
| Power Consumption- Water Heating | 5/15 KW |
| Power Consumption-Drying | 3.6 KW |
| Power Supply | AC220 V, 50 Hz , 380 V ,50Hz (Standard) 110V , 60 Hz (Optional) |
| Cleaning Rack | 3 pcs |
| Standard Cleaning Rack | BKIR36 x 3 |
| Standard Accessories | USB interface, One Bucket RBS A 156 leaner (1L), One Bucket RBS A 375 neutralizer (1L) |
| Optical Accessories | Cleaning Racks |
| Washer Chamber Size (WxDxH) | 600x629x820 mm |
| External Size (WxDxH) | 690x790x1970 mm |
| Package Size (WxDxH) | 790x940x2100 mm |
| Gross Weight | 320 kg |



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@biozef.com | Website: biozef.com