

# ELECTRIC BED 5 FUNCTION HBD11-05



# ELECTRIC BED 5 FUNCTION HBD11-05

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

## HBD11-05 ELECTRIC BED 5 FUNCTION

ABS head & foot boards can be removable

Aluminium alloy side rail, 6-piece, collapsed type with safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

Covered castors dia. 125mm

Powerful actuator system

With telescopic SS IV pole and previsions



## SPECIFICATIONS

Model	HBD11-05
Function	5
Back-rest lifting angle	80°±5°
Foot-rest lifting angle	40°±5°
Trendelenburg	12°
Reverse trendelenburg	12°
Height adjustment	500-750 mm (without mattress)
Drainage bag holders	Yes
IV pole	Yes
Dimension	2080x970x500 mm



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)