

# ELECTRIC BED 3/5 FUNCTION HBD13-05



# ELECTRIC BED 3/5 FUNCTION HBD13-05

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

## HBD13-05 ELECTRIC BED 3/5 FUNCTION

ABS head & foot boards can be removable

ABS side rail, 4-piece, tuck away type with safe lock

Bed sheet top made from perforated steel, epoxy coated and baked finish, with soft connection

Bed framework made from mild steel ERW rectangular tubes

Central locking foot pedal control with casters dia 150 mm

Powerful actuator system

With telescopic SS IV pole and previsions



## SPECIFICATIONS

Model	HBD13-05
Function	5
Back-rest lifting angle	$80^{\circ} \pm 5^{\circ}$
Foot-rest lifting angle	$40^{\circ} \pm 5^{\circ}$
Trendelenburg	$12^{\circ}$
Reverse trendelenburg	$12^{\circ}$
Height adjustment	500-750 mm (without mattress)
Drainage bag holders	Yes
IV pole	Yes



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)