

# THREE FUNCTION ELECTRIC BED HBD14-03



# THREE FUNCTION ELECTRIC BED HBD14-03

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

## HBD14-03 THREE FUNCTION ELECTRIC BED

Steel Framework, top class ABS & PP plastic

3 functions controlled by remote control

Castor with central brake system, wheel diameter: 150mm



### SPECIFICATIONS

Model	HBD14-03
Function	3
Back-rest lifting angle	0°-70°±5°
Foot-rest lifting angle	0°-30°±5°
Load Bearing	< 250 kg
Patient surface	1935Lx900W mm
Height adjustment	440-655 mm
IV pole	Yes
Dimension	2160Lx1120Wx440-655 mm
Net Weight	139 kg



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)