

# FIVE FUNCTION ELECTRIC BED HBD15-05



# FIVE FUNCTION ELECTRIC BED HBD15-05

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

## HBD15-05 FIVE FUNCTION ELECTRIC BED

Steel Framework, top class ABS & PP plastic

Functions activated by remote controller, with backup battery

Protective railing with angle indicators

Castor with central brake system, wheel diameter: 150mm



### SPECIFICATIONS

Model	HBD15-05
Function	Multi
Back-rest lifting angle	0°-65°±5°
Foot-rest lifting angle	0°-40°±5°
Trendelenburg	12°
Reverse trendelenburg	12°
Load Bearing	< 250 kg
Patient surface	1945Lx900W mm
Height adjustment	485-755 mm
IV pole	Yes
Dimension	2100Lx1050Wx485-755 mm
Net Weight	139 kg



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)