

MULTI- FUNCTION ELECTRIC BED HBD16-M



MULTI- FUNCTION ELECTRIC BED HBD16-M

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

HBD16-M MULTI- FUNCTION ELECTRIC BED

Steel Framework, top class ABS & PP plastic

Functions activated by electric controlling adopt luxurious membrane switch, easily operated and pressed comfortably

With the "reset" button to make the bed be flat from any position

A lock button with indicator light to show work state

Protective railing with angle indicators

castor with central brake system, wheel diameter: 150mm



SPECIFICATIONS

Model	HBD16-M
Function	Multi
Back-rest lifting angle	0°-65°±5°
Foot-rest lifting angle	0°-40°±5°
Trendelenburg	12°
Reverse trendelenburg	12°
Load Bearing	< 250 kg
Patient surface	1945Lx900W mm
Height adjustment	485-755 mm
IV pole	Yes
Dimension	2100Lx1050Wx485-755 mm
Net Weight	139 kg



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com