

# **ELECTRIC 3 FUNCTION MEDICAL BED HBD1B-03**



# ELECTRIC 3 FUNCTION MEDICAL BED HBD1B-03

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

## HBD1B-03 ELECTRIC 3 FUNCTION MEDICAL BED

ABS head & foot boards can be removable. With bumper wheels, 12 kg/pair, color option: blue, green with patient name card holder

2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

Central locking S.S foot pedal with casters dia. 5"

With 3 motors from Taiwan. With battery.

Backrest & knee-rest up/down at same time

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes

Optional Accessories: Manual CPR/electric CPR



## SPECIFICATIONS

Model	HBD1B-03
Function	3
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Height adjustment	430-720 mm (without mattress)
IV pole	Yes
Dimension	2160Lx980Wx430-720 mm



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)