

ELECTRIC 2 FUNCTION MEDICAL BED HBD1C-02



ELECTRIC 2 FUNCTION MEDICAL BED HBD1C-02

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

HBD1C-02 ELECTRIC 2 FUNCTION MEDICAL BED

ABS head & foot boards can be removable. with bumper wheels, 12kg/pair, color option: brown, blue, green, with patient name card holder

ABS side rail, 4-piece, tuck away type with safe lock. with angle display

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

ABS foot pedal with casters dia. 5"

With 2 Linak motors from Denmark. with battery

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes



SPECIFICATIONS

Model	HBD1C-02
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
IV pole	Yes
Dimension	2260Lx1060Wx500H mm



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com