

# **ELECTRIC 5 FUNCTION MEDICAL BED HBD1D-05**



# ELECTRIC 5 FUNCTION MEDICAL BED HBD1D-05

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

## HBD1D-05 ELECTRIC 5 FUNCTION MEDICAL BED

ABS head & foot boards can be removable. 6.3 kg/pair, color options: blue with patient name card holder

ABS side rail, 4pcs, tuck away type with safe lock-out

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

Central locking S.S foot pedal with casters dia. 5"

With 4 motors from Taiwan

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes



### SPECIFICATIONS

Model	HBD1D-05
Function	5
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Trendelenburg	0°~12°
Reverse trendelenburg	0°~12°
Height adjustment	430-720 mm (without mattress)
IV pole	Yes
Dimension	2120Lx1030Wx430-720 mm



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)