

ELECTRIC 3 FUNCTION MEDICAL BED HBD1F-03



ELECTRIC 3 FUNCTION MEDICAL BED HBD1F-03

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

HBD1F-03 ELECTRIC 3 FUNCTION MEDICAL BED

ABS head & foot boards can be removable. With bumper wheels, 12 kg/pair, color option: brown, blue, green with patient name card holder

ABS side rail, 4pcs, tuck away type with safe lock. With angle display

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

Central locking both side foot pedal with casters dia. 5"

With 3 linak motors from Denmark. With battery.

Backrest & knee-rest up/down at same time

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes

Optional Accessories: CPR/electric CPR



SPECIFICATIONS

Model	HBD1F-03
Function	3
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Height adjustment	430-720 mm (without mattress)
IV pole	Yes
Dimension	2260Lx1060Wx430-720 mm



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com