

# 5 FUNCTIONS ELECTRIC IUC BED HBD1G-05



# 5 FUNCTIONS ELECTRIC IUC BED HBD1G-05

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

## HBD1G-05 5 FUNCTIONS ELECTRIC IUC BED

With one-key electric CPR release functions

With manual CPR functions

With one-key TR functions

With one-key examination position

With one-key chair position

Frame: made of cold-rolled steel plate, treated by electro-coating & powder coating

4 pcs detachable PP side rails

Build-in control embedded inside the railing panel

Quiet and robust electric actuators control of all reliable operation

Nursing screen remote handset control of all electrically operated functions and movements

Pedal central-controlled caster brake system, dia. 125 mm

Advanced medical motor system (4pcs Linak motors, 1 pcs control box, 1pcs handset)

Accessories: 1pc telescopic SS I.V. pole; 4 pcs urine bag hooks; 4pcs I.V. holes



## SPECIFICATIONS

|                         |                               |
|-------------------------|-------------------------------|
| Model                   | HBD1G-05                      |
| Function                | 5                             |
| Back-rest lifting angle | 0°-80°±5°                     |
| Foot-rest lifting angle | 0°-40°±5°                     |
| Trendelenburg           | 0°~12°                        |
| Reverse trendelenburg   | 0°~12°                        |
| Height adjustment       | 430-720 mm (without mattress) |
| IV pole                 | Yes                           |
| Dimension               | 2260Lx1060Wx430-720 mm        |
| Net Weight              | 250 kg                        |



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)