

4 CRANK/ MANUAL 5 FUNCTION MEDICAL BED HBD1J-05



4 CRANK/ MANUAL 5 FUNCTION MEDICAL BED

HBD1J-05

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

HBD1J-05 4 CRANK/ MANUAL 5 FUNCTION MEDICAL BED

ABS head & foot boards can be removable. with bumper, 8kg/pair, color option: brown, blue, with patient name card holder

2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

With casters dia. 5" with individual brake

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes

Optional Accessory: Shoe shelf, mosquito net rack, dinner table



SPECIFICATIONS

| | |
|-------------------------|-------------------------------|
| Model | HBD1J-05 |
| Function | 5 |
| Back-rest lifting angle | 0°-80°±5° |
| Foot-rest lifting angle | 0°-40°±5° |
| Trendelenburg | 0°~12° |
| Reverse trendelenburg | 0°~12° |
| Height adjustment | 500~780 mm (without mattress) |
| IV pole | Yes |
| Dimension | 2150Lx1060Wx500-780 mm |



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com