

AUTOMATIC BED HBD1K-A



AUTOMATIC BED HBD1K-A

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

HBD1K-A AUTOMATIC BED

4 pcs of European style guardrail, upside down structure

Luxury central control castors, brake the bed with one stamp

With 4 bumpers



SPECIFICATIONS

Model	HBD1K-A
Back-rest lifting angle	0°-65°±5°
Foot-rest lifting angle	0°-35°±5°
Trendelenburg	12°
Reverse trendelenburg	12°
Dimension	2200x1020x420-780 mm



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com