

AUTOMATIC BED HBD1K-A



AUTOMATIC BED HBD1K-A

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

HBD1K-A AUTOMATIC BED

4 pcs of European style guardrail, upside down structure
Luxury central control castors, brake the bed with one stamp
With 4 bumpers



SPECIFICATIONS

| | |
|-------------------------|----------------------|
| Model | HBD1K-A |
| Back-rest lifting angle | 0°-65°±5° |
| Foot-rest lifting angle | 0°-35°±5° |
| Trendelenburg | 12° |
| Reverse trendelenburg | 12° |
| Dimension | 2200x1020x420-780 mm |



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com