

AUTOMATIC BED HBD1L-A1



AUTOMATIC BED HBD1L-A1

An electric Hospital Bed is ideal if you're searching for a bed that can be completely controlled with the click of a button. Features with integrated multi-function and can fulfil the needs of various operations. Electric hospital bed allows height adjustment, as well as providing optimal height for caregivers to work with patients; raising of the head, including into Fowler's position, to allow for feeding, ease of breathing and other benefits resulting from reclining upright; and raising of the feet to help with ease of movement and blood flow.

Used in Medical unit, ICU.

Also known as Medical bed, hospital bed, Hospital cot.

HBD1L-A1 AUTOMATIC BED

With the nurse controller and CPR function

Bed board adopts HPL board, x-ray photograph for whole body is available

2 step adjustment at shank section



SPECIFICATIONS

Model	HBD1L-A1
Back-rest lifting angle	0°-65°±5°
Foot-rest lifting angle	0°-35°±5°
Bed panel extention	Manual/180mm
Load Bearing	240 kg
Patient surface	1950x910 mm
Dimension	2170x1080x400~750 mm



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com