

MANUAL BED DOUBLE CRANK HBD24-MC



MANUAL BED DOUBLE CRANK HBD24-MC

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance

Used in Recuperating, Nursing.

Also known as Medical bed, hospital bed, Hospital cot, Electric ICU bed.

HBD24-MC MANUAL BED DOUBLE CRANK

ABS head and boards can be removable

ABS side rail, 4-piece, tuck away type with safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

Central locking foot pedal control with castors dia 150 mm

With shoes holder

With telescopic SS IV pole and previsions



SPECIFICATIONS

Model	HBD24-MC
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
IV pole	Yes
Dimension	2080x970x500 mm



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com