

# MANUAL BED 3 CRANK HBD26-03



# MANUAL BED 3 CRANK HBD26-03

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance

Used in Recuperating, Nursing.

Also known as Medical bed, hospital bed, Hospital cot, Electric ICU bed.

## HBD26-03 MANUAL BED 3 CRANK

ABS head and boards can be removable

ABS side rail, 4-piece, tuck away type With safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish, with soft connection.

Bed framework made from mild steel ERW rectangular tubes.

Central locking foot pedal control with castors dia 150 mm.

With drainage bag holders

With telescopic SS IV pole and previsions



## SPECIFICATIONS

Model	HBD26-03
Function	3
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Height adjustment	500-750 mm (without mattress)
Drainage bag holders	Yes
IV pole	Yes
Dimension	2180x1030x500 mm



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)