

# TWO CRANK MANUAL BED HBD28-02



# TWO CRANK MANUAL BED HBD28-02

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance

Used in Recuperating, Nursing.

Also known as Medical bed, hospital bed, Hospital cot, Electric ICU bed.

## HBD28-02 TWO CRANK MANUAL BED

Steel framework, top class ABS & PP plastic

Collapsible Aluminium alloy protective railing; Detachable PP head & foot board

4 pcs universal castor, dia. 125 mm, Individual brake

A shoe rack under the bed board



## SPECIFICATIONS

Model	HBD28-02
Function	2
Back-rest lifting angle	0°-75°±5°
Foot-rest lifting angle	0°-45°±5°
Patient surface	1910Lx830W mm
Height adjustment	470 mm
IV pole	Yes
Dimension	2080Lx980Wx470 mm
Net Weight	65 kg



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)