

3 CRANK MANUAL BED HBD29-03



3 CRANK MANUAL BED HBD29-03

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance

Used in Recuperating, Nursing.

Also known as Medical bed, hospital bed, Hospital cot, Electric ICU bed.

HBD29-03 3 CRANK MANUAL BED

Steel framework, top class ABS & PP plastic

PP head & foot board and bed board, detachable and easy to clean 6 columns
collapsible Aluminium alloy protective railing

4 pcs universal castor, dia. 125 mm, Individual brake



SPECIFICATIONS

Model	HBD29-03
Function	3
Back-rest lifting angle	0°-65°±5°
Foot-rest lifting angle	0°-45°±5°
Patient surface	1930x830 mm
Height adjustment	750-575 mm
IV pole	Yes
Dimension	2160x960x530- 730 mm



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com