

# 4 CRANK MANUAL BED HBD2A-04



# 4 CRANK MANUAL BED HBD2A-04

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance

Used in Recuperating, Nursing.

Also known as Medical bed, hospital bed, Hospital cot, Electric ICU bed.

## HBD2A-04 4 CRANK MANUAL BED

Steel framework, top class ABS & PP plastic

PP material head / foot board, detachable and easy-clean

4 pcs drop-down PP protective railing, controlled by gas spring

Central brake castors, dia. 150 mm. Brake pedal locate below the foot board



## SPECIFICATIONS

Model	HBD2A-04
Function	5
Back-rest lifting angle	0°-60°±5°
Foot-rest lifting angle	0°-40°±5°
Trendelenburg	0-12°
Reverse trendelenburg	0-12°
Patient surface	1950x900 mm
Height adjustment	755-485 mm
IV pole	Yes
Dimension	2200Lx1050Wx485-755 mm
Net Weight	117 kg



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)