

# MANUAL 2

# FUNCTION MEDICAL

# BED HBD2J-02



# MANUAL 2 FUNCTION MEDICAL BED HBD2J-02

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance

Used in Recuperating, Nursing.

Also known as Medical bed, hospital bed, Hospital cot, Electric ICU bed.

## HBD2J-02 MANUAL 2 FUNCTION MEDICAL BED

ABS head & foot boards can be removable. with bumper wheels, 12kg/pair, color option: brown, blue, green, with patient name card holder

ABS side rail, 4-piece, tuck away type with safe lock. with angle display

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes.

Central locking both side foot pedal with casters dia. 5"

Manual 2 functions adjusted by hidden type hand crank

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes



## SPECIFICATIONS

Model	HBD2J-02
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
IV pole	Yes
Dimension	2220Lx1060Wx500 mm



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)