

# SINGLE CRANK HOSPITAL BED HBD2M-HM



# SINGLE CRANK HOSPITAL BED HBD2M-HM

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance

Used in Recuperating, Nursing.

Also known as Medical bed, hospital bed, Hospital cot, Electric ICU bed.

## HBD2M-HM SINGLE CRANK HOSPITAL BED



### SPECIFICATIONS

Model	HBD2M-HM
Load Bearing	240 kg
Patient surface	2150x950 mm
Height adjustment	500 mm (from floor)



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)