

MANUAL HOSPITAL BED WITH 2 FUNCTION HBD2P- HM2



MANUAL HOSPITAL BED WITH 2 FUNCTION

HBD2P-HM2

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance

Used in Recuperating, Nursing.

Also known as Medical bed, hospital bed, Hospital cot, Electric ICU bed.

HBD2P-HM2 MANUAL HOSPITAL BED WITH 2 FUNCTION



SPECIFICATIONS

Model	HBD2P-HM2
Back-rest lifting angle	0°-75°
Foot-rest lifting angle	0°-35°
Load Bearing	240 kg
Patient surface	1950x830 mm
Dimension	2200x1020x500 mm



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com