

DELIVERY BED HBD32-DB1



DELIVERY BED HBD32-DB1

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.

Used in Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation.

Also known as Gyneacological table, Gyneacological bed, Examination table, Examination bed, Delivery bed, Gynaecological Delivery Bed.

HBD32-DB1 DELIVERY BED

Table frame made of tubular stainless steel; Fixed bases with protective plastic feet

Adjustable three sections with independent controls and trendelenburg positioning

Bed top is high-density polyurethane foam, and covered with washable faux leather: waterproof, flame retardant

Mattress thickness: 40 mm

One pair of padded shoulder support, one pair of adjustable knee support; one pair of handgrip

With slide-out S.S bowl

Knock down construction



SPECIFICATIONS

Model	HBD32-DB1
Load Bearing	150 kg
Dimension	1800x600x800 mm



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com