

# MANUAL OBSTERIC BED HBD33-MO



# MANUAL OBSTERIC BED HBD33-MO

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.

Used in Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation. Also known as Gyneacological table, Gyneacological bed, Examination table, Examination bed, Delivery bed, Gynaecological Delivery Bed.

## HBD33-MO MANUAL OBSTERIC BED

The main bed and auxiliary bed are connected with each other, four foot wheels are equipped on auxiliary bed and brake device is spared

Foot step hydraulic elevating, Imported Y type sealing ring and is durable

The mattress adopts high density memory sponge and is once formed, easily cleaned and sterilized

Rotary brake device which is easy for moving operating table

Base seat and upright column covering are made of high quality 304 stainless steel



## SPECIFICATIONS

Model	HBD33-MO
Back-rest lifting angle	$\geq 50^\circ$
Back-rest downwards	$\geq 10^\circ$
Seat board upwards	$\geq 30^\circ$
Seat board downwards	$\geq 10^\circ$
Working area	650x600
Dimension	1950Lx600Wx750~950 mm



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)